

DIONISOS:

A new experiment studying the
dynamics of plasma-surface interactions

Dennis Whyte, Graham Wright, Alan Hoskinson
University of Wisconsin - Madison

PFC Meeting

May 2004

University of Illinois

Our present lack of confidence in PSI issues primarily arises from the poor diagnosis of PSI effects

- Most of the controlling PSI physics is in hand or in a “mature” research stage
 - Physical sputtering, chemical sputtering
 - Plasma edge sheath theory
 - Atomic physics (ionization, etc.)
 - *Notable exception: \perp particle transport in edge plasmas.*
- By the same token, edge plasma diagnosis is mature.
- However we have essentially no direct diagnosis of how & when plasma-facing surfaces are being modified by the plasma.

The PSI diagnostic challenge

- The plasma edge is dominated by its own self-consistent recycling pattern of fuel and impurities
 - Easily implemented plasma diagnosis cannot inform us about the net surface effects (H retention, erosion, etc.)
- The lack of direct surface diagnosis cripples our ability to understand the net effects of PSI.
 - We have largely relied on surface “archeology” to inform us about surface modifications. This is unacceptable; the equivalent of basing ITER burn predictions on one-year integrated neutron fluence from JET!!)
 - Dedicated ex-situ PSI experiments (e.g. DiMES) very expensive to implement (run-time) and covers $\ll 1\%$ of wall surfaces.
- The challenge is to develop reliable in-situ surface diagnostic techniques that can operate in real-time and be placed in many locations.

In-situ PSI diagnosis topics

1. Erosion & deposition using quartz crystal microbalances
2. DIONISOS: Accelerator-based ion beam surface analysis experiment being developed at UW-Madison
3. Radio-isotope alpha emission for remote, in-situ “ion beam analysis” in confinement experiments.

QMB: Quartz crystal microbalances

- Measures net mass change on surface coating of a crystal by changed resonance frequency.
- **Advantages:**
 - Extremely sensitive: can measure ~mono-layer net changes in surface
 - Large dynamic range: 0.1-> 1000 nm
 - Can measure both deposition and erosion of a pre-deposited film.
 - Commercially available technology.
- **Drawbacks:**
 - No element discrimination.
 - Cannot tolerate any significant heat flux: small thermal mass.
 - Highly temperature / environment sensitive
 - Limited absolute range: ~ 1 micron << expected deposit layers.


QMB: Quartz crystal microbalances

- **Implementation needs:**
 - Special design for in-vessel electronics (ceramic based) to tolerate high baking temperatures (~ 300 C)
- **Present status.**
 - QMB successfully implemented in JET divertor...clearly demonstrating need for real-time PSI diagnosis.
 - UW-Madison, MIT and GA now collaborating on implementation of QMBs in Alcator C-Mod (Mo/B) and DIII-D (C) tokamaks.
- **Summary:** Primarily due to their intolerance to heat flux, QMBs are relegated to measuring erosion/redeposition in “hidden” areas (baffles, main-wall ports, behind tiles, etc.). Going into new DIII-D lower divertor in tile gaps.

Ion Beam Analysis of Surfaces


- Interaction of \sim MeV ions with surface atoms is accurately described from well-known physics
 - Rutherford elastic scattering ($M_{\text{target}} > M_{\text{projectile}}$)
 - Forward elastic recoil ($M_{\text{target}} < M_{\text{projectile}}$)
 - Inelastic nuclear reactions (all reaction σ known)
- Exploit this understanding to diagnose near-surface (\sim 1-10 microns) properties of materials.
- **Advantages:**
 - Great flexibility: element/isotope sensitivity, erosion/deposition/H retention available in single diagnostic.
 - Measurement accuracy over appropriate depth range for PSI.
 - Proven application in fusion PSI (see examples from DiMES)
 - Simple solid-state detectors for scattered charged particles
- **Drawbacks:**
 - Need \sim MeV energy ions

DiMES examples of Ion Beam Analysis of Surfaces




DiMES examples of Ion Beam Analysis of Surfaces

- Carbon erosion/deposition is determined by the change in depth of an implanted silicon marker measured by 2 MeV helium Rutherford backscattering (RBS). Detection limit $\sim \pm 10$ nm.
- Metal erosion is determined from the change in thickness of thin metal films measured by RBS (W,V).
- Nuclear reaction yields give areal density of:
Deuterium $D(^3\text{He},p)\alpha$
Boron $^{11}\text{B}(p,\alpha)^8\text{Be}$





Near-surface transport

H Retention



Low-Z erosion



DIONISOS:

Dynamics of IONic Implantation & Sputtering On Surfaces

- A new experimental facility being constructed at UW-Madison.
- *Goal:* Accurately measure the real-time response of the PFC material surface to plasma bombardment using in-situ high-energy ion beam analysis.
- Design philosophy:
 - Exploit ion beam analysis (IBA) techniques
 - Mature analysis tool widely used in R&D.
 - Vastly different ion energies (30 eV vs. > 1 MeV) allow for simultaneous surface interaction with little interference.
 - Use previously developed axisymmetric plasma sources
 - *Helicon:* Steady-state plasma /w good density control.
 - *Plasma gun:* high density/flux pulsed plasma source.
 - Initially focus on fusion experiments, but generic PSI tool.



**Principal component of DIONISOS is now operational at
UW-Madison: *1.7 MV tandem ion accelerator***



MODEL 5SDH-4 PELLETRON ACCELERATOR


- Features:
 - Dual sources (sputtering and RF plasma), > 100 beam species available.
 - High energy (≤ 10 MeV for higher Z beams)
 - High current beamlines (> 0.1 mA) for implantation and irradiation.

DIONISOS Experiment Setup




DIONISOS Experiment Setup

**MeV Ion beam
In-situ implantation &
surface diagnosis**



Experimental R&D: Development of plasma compatible ion beam analysis detector assemblies



- Serves purpose of protecting s.s. detector from
 - Plasma light emission.
 - Low energy plasma ions.
 - Sputtered and CX neutrals.

Three proposed areas of study for DIONISOS


- **Measurement and modeling of near-surface cross-field ion transport**
 - **Determine controlling parameters for the magnitude and locations of net erosion / redeposition of PFC relevant to fusion (C, W, Mo)**
- **The dynamic release of fuel and impurity particles from surfaces exposed to transient, high-density plasmas.**
- **The dynamics of hydrogenic / tritium fuel trapping in plasma-deposited films, for single and multiple species materials.**

Plasma sources for DIONISOS


- **Helicon plasma source (collaboration with N. Hershkowitz)**
 - **Steady-state ionizing plasma**
 - **Solenoid field convenient linear control of plasma density**
 - **H, D, Ne, Ar**
- **Plasma gun (collaboration with G. Fiksel, C. Forest)**
 - **Developed for helicity injection on MST (Fiksel et al).**
 - **Pulsed ($\sim 1\text{-}10$ ms) H/D plasmas with ~ 60 s rep rate.**
 - **$\sim 100\%$ ionization, $T_e \sim 20$ eV high density ($\sim 10^{20} \text{ m}^{-3}$)**
 - **$\sim 100 \text{ MW/m}^2$ ($10 \text{ MJ/m}^2/\text{s}^{1/2}$) per 100 V bias approaches ablation/melt limits.**
 - **Capable of current densities $\sim \text{kA} / \text{cm}^2$**

Ion beam surface analysis tools are operational on DIONISOS

- **Rutherford Backscattering Spectroscopy**
- **Elastic Recoil Detection (H detection)**
- **Nuclear Reaction Analysis.**
- **Particle Induced Gamma Emission.**
 - Carbon-13 detection with ~2 nm depth resolution for DIII-D experiment.




Schematic for IBA diagnosis of high-Z PFC materials




- RBS measures net erosion rate of high-Z film to ~ 5 nm.
- Forward recoil (ERD) provides depth-resolved H/D concentration.
- Highly resolved spatial profiles of isotope resolved H species with NRA.

Schematic for IBA diagnosis of low-Z PFC materials



- RBS measures net erosion rate of bulk material by changes in high-Z marker previously implanted with ion beam.
- ERD and NRA provide real-time hydrogenic concentrations and diffusion in deposited films

Schematic for IBA diagnosis of H dynamics in PFC




- Exploit large forward recoil cross-section for high mass projectiles
- Same ion-PFC species allows for high current, fast hydrogenic diagnosis
 - Case shown: ERD using C-H recoil can measure H profile in ~ 10 ms time
- Allows for dynamic study of H implantation, diffusion and release under transient conditions.

Remote IBA using radio-isotope alpha emitters

- While IBA is the tool of choice for surface diagnosis, how can we provide real-time measurements in a confinement device where we have no access with ion beam?
- Produce an ~mono-energetic alpha ion “beam” using natural alpha emitters
 - Present focus on Po-210 (138 day half-life, $E_{\alpha}=5.4$ MeV)
 - Balances need for alpha flux vs. diagnostic lifetime.




Schematic of Alpha Surface Analyzers (ASA) design for real-time PSI diagnosis in confinement devices



- **Exploit intrinsic magnetic field to rotate surface in/out of contact with main PFC surface.**
 - Large thermal mass and contact allow for surfaces in high heat flux areas.
 - Controls exposure duration of surface in each discharge.
- **In analysis position, the alpha scattering is detected by s.s. detectors for**
 - RBS: net erosion/deposition and surface stoichiometry
 - ERD: Hydrogenic retention over ~ 5 microns.
- **Po-210 source will allow for surface diagnosis between each shot (~10 minute acquisition time).**

First tests of Alpha Surface Analyzers (ASA) design are being carried out for Alcator C-Mod



- Multi-layer B/Mo films allow for accurate measurements of net erosion, net deposition and mixed-material analysis in C-Mod.
- Other R&D
 - Po-210 radiochemistry for film deposition.
 - Detector geometry optimization.